

SPECULUMS, RETRACTORS AND FIXATION RINGS

Wire Speculums.....	1 – 1-3
Adjustable Speculums.....	1 – 4
Temporal Speculums Lieberman Style	1 – 5
Nasal Speculums Lieberman Style.....	1 – 6
Solid Blade Speculums Lieberman Style.....	1 – 7
Triple Post Speculums	1 – 8
Iris Retractors	1 – 9
Lachrymal Sac Retractors.....	1 – 9-10
Lid Retractors	1 – 10
Orbital Retractors.....	1 – 11
Tissue Retactors	1 – 11
Flieringa Scleral Fixation Rings	1 – 12
Fixation Rings	1 – 12-13
Fixation Hooks	1 – 14

SPECULUMS

WIRE SPECULUMS

Barraquer Wire Speculums



AS1-001 HW


- Adult Size
- 14mm Blades
- Heavy Wire
- Stainless Steel

AS1-001S

- Adult Size
- 14mm Blades
- Standard Wire
- Stainless Steel

AS1-001T

- Adult Size
- 14mm Blades
- Standard Wire
- Titanium



AS1-002 HW

- Child Size
- 11mm Blades
- Heavy Wire
- Stainless Steel

AS1-002S

- Child Size
- 11mm Blades
- Standard Wire
- Stainless Steel

AS1-002T

- Child Size
- 11mm Blades
- Standard Wire
- Titanium



AS1-004S

- Neonatal Size
- 4mm Blades
- Standard Wire
- Stainless Steel

AS1-004T

- Neonatal Size
- 4mm Blades
- Standard Wire
- Titanium

Barraquer Speculums With Solid Blades



AS1-005S

- Adult Size
- 14mm Blades
- Standard Wire
- Stainless Steel

AS1-005T

- Adult Size
- 14mm Blades
- Standard Wire
- Titanium

AS1-006S

- Child Size
- 11mm Blades
- Standard Wire
- Stainless Steel

AS1-006T

- Child SizeT
- 11mm Blades
- Standard Wire
- Titanium



AS1-007S


- Neonatal Size
- 4mm Blades
- Standard Wire
- Stainless Steel

AS1-007T

- Neonatal Size
- 4mm Blades
- Standard Wire
- Titanium

SPECULUMS

Kratz-Barraquer Wire Speculums



AS1-008

- Adult Size
- Open Blades
- 14mm Blades
- Standard Wire

AS1-009

- Adult Size
- Open Blades
- 14mm Blades
- Heavy Wire

AS1-010

- Child Size
- Open Blades
- 11mm Blades
- Standard Wire

Nasal Barraquer Wire Speculums



AS1-011

- Adult Size
- 14mm Blades
- Standard Wire

AS1-012

- Adult Size
- 14mm Blades
- Heavy Wire

AS1-013

- Child Size
- 11mm Blades
- Standard Wire

AS1-014

- Neonatal Size
- 4mm Blades
- Standard Wire

Nasal Barraquer Speculums



AS1-015

- Adult Size
- Solid Blades
- 14mm Blades

AS1-016

- Child Size
- Solid Blades
- 11mm Blades
- Standard Wire

AS1-017

- Neonatal Size
- Solid Blades
- 4mm Blades
- Standard Wire

Nasal Kratz-Barraquer Wire Speculums



AS1-018

- Adult Size
- Open Blades
- 14mm Blades
- Standard Wire

AS1-019

- Adult Size
- Open Blades
- 14mm Blades
- Heavy Wire

AS1-022

- Child Size
- Open Blades
- 11mm Blades
- Standard Wire

SPECULUMS

Alphonso Eye Speculums



AS1-023S

- With Discs
- Newborn Size
- 5mm and 8mm Blades
- 27mm Blade Spread
- 45mm Overall Length
- Stainless Steel

AS1-023T

- With Discs
- Newborn Size
- 5mm and 8mm Blades
- 27mm Blade Spread
- 45mm Overall Length
- Titanium

Jaffe Lid Retractors

AS1-024

- Pair
- 15mm
- Adult

Mouse Speculum



AS1-040

- 4mm Blade
- 22mm Overall Length

Neyvas Drape Speculum

AS1-080

- Open Wire
- 16mm Blades
- Designed to retain surgical drape under the eyelids

SPECULUMS

ADJUSTABLE SPECULUMS

Cook Eye Speculums



AS1-050

- Small
- With Locking Screw
- Stainless Steel

AS1-051

- Medium
- With Locking Screw
- Stainless Steel

AS1-052

- Large
- With Locking Screw
- Stainless Steel

Sauer Lid Speculums



AS1-053

- Child
- 8mm Blades
- 20mm Blade Spread
- 35mm Overall Length
- Titanium

AS1-054

- Newborn
- 5mm Blades
- 19mm Blade Spread
- 26mm Overall Length
- Titanium

AS1-055

- Premature
- 4mm Blades
- 18mm Blade Spread
- 24mm Overall Length
- Titanium

Frankel Infant Lid Speculum



AS1-060

- 10mm Blades
- Adjustable
- Stainless Steel

Ong Glaucoma Lid Speculum

AS1-115LT

- For use in the LEFT Eye
- Larger Inferior Blade pushes on the Inferior Conjunctival Fornix to turn eye down.

AS1-115RT

- For use in the RIGHT Eye
- Larger Inferior Blade pushes on the inferior Conjunctival Fornix to turn eye down.

SPECULUMS

TEMPORAL SPECULUM LIEBERMAN STYLE

Lieberman Temporal Speculum



AS1-100VS

- Adult Size
- V Blades
- 14mm Blades
- Stainless Steel

AS1-100VT

- Adult Size
- V Blades
- 14mm Blades
- Titanium

AS1-100CVT

- Child Size
- V Blades
- 11mm Blades
- Titanium

AS1-100CVT-COVAR

- Especially Designed for Small Eyes
- V Blades
- 11mm Blades
- Titanium
- Overall Size 65mm

AS1-100CVAS

- Child Size
- V Blades
- Stainless
- Aspirating



AS1-100KS

- Adult Size
- Kratz Blades
- 14mm Blades
- Stainless Steel

AS1-100KT

- Adult Size
- Kratz Blades
- 14mm Blades
- Titanium



AS1-100RT

- Adult Size
- Round Blades
- 14mm Blades
- Especially Designed for LASIK
- Titanium

SPECULUMS

NASAL SPECULUM LIEBERMAN STYLE

Lieberman Nasal Speculums



AS1-110VT

- Adult Size
- V Blades
- 14mm Blades
- Titanium

AS1-110CVT

- Child Size
- V Blades
- 11mm Blades
- Titanium



AS1-110KT

- Adult Size
- Kratz Style Blades
- 14mm Blades
- Titanium



AS1-110RT

- Round Blades
- Especially Designed for LASIK
- 14mm Blades
- Titanium

SPECULUMS

SOLID BLADE SPECULUM LIEBERMAN STYLE

Kershner Reversible Solid Blade Speculum




AS1-120S

- Solid Blades
- Temporal and Nasal Fit
- Stainless Steel

AS1-120T

- Solid Blades
- Temporal and Nasal Fit
- Titanium

Kershner Reversible Speculum with Aspiration



AS1-120AS

- Solid Blades
- Supplied with Silicone Tubing and Adaptor
- Temporal and Nasal Fit
- Stainless Steel

AS1-120AT

- Solid Blades
- Supplied with Silicone Tubing and Adaptor
- Temporal and Nasal Fit
- Titanium

Kershner Reversible Speculum Open Blades



AS1-121T

- 20mm Closed Wire Blades
- Titanium

SPECULUMS

TRIPLE POST SPECULUMS

Triple Post Adjustable Speculums



AS1-130S

- Stainless Steel

AS1-130T

- Titanium

FEMTO SPECULUMS

Aurora Femto Speculums



AS1-140

- 15mm Adjustable "K" Shaped Blades
- 72mm Overall Length
- Titanium

Lieberman Femto Speculum



AS1-141KS

- Lieberman Femtosecond Adjustable Lid Speculum
- Wire Silicone Protected Blades
- Temporal Model
- 14mm Open Square Blades with Flares
- Stainless Steel

AS1-141VS

- Lieberman Femtosecond Adjustable Lid Speculum
- Wire Silicone Protected Blades
- Temporal Model
- 14mm Open V Blades with Flares
- Stainless Steel

Crozafon Femto Speculum



AS1-142

- Crozafon Femtosecond Adjustable Slide Speculum with Lock Control Screw
- Angled Shafts
- 20mm Open Round Blades with Flares
- Stainless Steel

RETRACTORS

IRIS RETRACTORS & CAPSULE SUPPORT SYSTEM

Disposable Iris Retractors



AS1-203

- 1 Pack of 4 pcs (4 retractors)

AS1-205

- 1 Pack of 5 pcs (5 retractors)



Reusable Iris Retractor



AS1-207

- 1 Pack of 6 Retractors in Teflon Container



McKool Cataract Support System

AS1-225-5

- Box of 5

LACHRYMAL SAC RETRACTORS

Stevenson Lachrymal Sac Retractor

AS1-230



RETRACTORS

Knapp Lachrymal Sac Retractor



AS1-231

- 2 Sharp Prongs



AS1-232

- 3 Blunt Prongs

AS1-233

- 3 Sharp Prongs



AS1-234

- 4 Blunt Prongs

AS1-235

- 4 Sharp Prongs

LID RETRACTORS

Desmarres Lid Retractor



AS1-240

- Size 0, 8mm

AS1-241

- Size 1, 11mm

AS1-242

- Size 2, 14mm

AS1-243

- Size 3, 17mm

Jaeger Lid Plate



AS1-244S

- 20mm and 23mm Wide Blades
- Matte Finish
- Stainless Steel

AS1-244T

- 20mm and 23mm Wide Blades
- Matte Finish
- Titanium

RETRACTORS

ORBITAL RETRACTORS

Schepens Orbital Retractor



AS1-250

- 13mm x 33mm Blade
- 4.5mm Wide Notch in Blade
- Titanium

Arruga Globe Retractor



AS1-251

Double Ended Orbital Globe Retractor-Elevator



AS1-252

- 10mm / 14mm

AS1-253

- 16mm / 25mm

TISSUE RETRACTORS

Helveston Tissue Retractor



AS1-254

- “Barbie”
- 7mm Wide

AS1-255

- “Big Barbie”
- 9mm Wide

AS1-256

- “Great Big Barbie”
- 11mm Wide Concave Blade

FIXATION RINGS

FLIERINGA SCLERAL FIXATION RINGS

Flieringa Scleral Fixation Ring

AS1-300-Set

- Full set of 9 (14mm-22mm)



AS1-300-14

- 14mm Diameter



AS1-300-15

- 15mm Diameter



AS1-300-16

- 16mm Diameter

AS1-300-17

- 17mm Diameter

AS1-300-18

- 18mm Diameter

AS1-300-19

- 19mm Diameter

AS1-300-20

- 20mm Diameter

AS1-300-21

- 21mm Diameter

AS1-300-22

- 22mm Diameter

FIXATION RINGS

Fine/Thornton Fixation Ring With Swivel



AS1-350S

- 13mm Diameter
- 11mm Opening
- Stainless Steel

AS1-350T

- 13mm Diameter
- 11mm Opening
- Titanium



AS1-351S

- 16mm Diameter
- 14mm Opening
- Stainless Steel

AS1-351T

- 16mm Diameter
- 14mm Opening
- Titanium

Gimbel Stabilization Ring - Offset



AS1-352S

- 13mm with 9mm Opening
- Blunt Teeth on Both Sides of Swivel
- Stainless Steel

AS1-352T

- 13mm with 9mm Opening
- Blunt Teeth on Both Sides of Swivel
- Titanium

LASIK Hinge Protector Fixation Ring



AS1-353T

- 11mm Inner Diameter
- Swivel
- Titanium

FIXATION RINGS

Lu-Mendez Fixation Ring



AS1-354

- 12mm Inner Diameter
- 18mm Outer Diameter
- Blunt Teeth
- Stainless Steel

Mendez Degree Gauge With Fixation



AS1-355

- Atraumatic Teeth for Fixation
- Calibrated Marks Every 10 Degrees from 0-180
- Internal 13mm / External 19mm
- Titanium

Fixation Ring Nichamin Style



AS1-356

- Nichamin Fixation Ring and Gauge
- Open 14mm Ring with 8mm Cut-out
- Fine teeth
- Gauge Marked Every 10 Degrees
- Two Broad Marks to Align with Steep Axis
- Round handle

Globe Fixation Ring-3/4 Open Ring With Fixation

AS1-357

- 3/4 Open Ring
- Stainless Steel

Globe Fixation Ring-Full Ring

AS1-358

- Full Ring
- Stainless Steel

FIXATION HOOKS

FIXATION HOOKS

Double Fixation Hooks



AS1-380

- Small
- 1.5mm Spread
- Sharp
- 5"mm (127mm)
Overall Length

AS1-381

- Large
- 2.0mm Spread
- Sharp
- 5"mm (127mm)
Overall Length

Twist Fixation Hooks



AS1-382L

- Left Twist

AS1-383L

- Left Twist

AS1-382R

- Right Twist
- 0.3mm

AS1-383R

- Right Twist
- 0.5mm