

**BLADE BREAKERS, BLADE HOLDERS  
BURRS, SPUDS AND TREPHINES**

Blade Breakers and Holders .....	4 – 1
Knives and Spuds .....	4 – 2-3
Alger Brushes .....	4 – 3
Pterygium Instruments .....	4 – 4
Corneal Trephines .....	4 – 5

## BLADE BREAKERS and HOLDERS

### BLADE BREAKERS AND HOLDERS

#### Barraquer Blade Breaker and Holder



**AS4-001S**

• Stainless Steel

**AS4-001T**

• Titanium

#### Castroviejo Blade Breaker and Holder



**AS4-002**

#### Bard Parker Handle



**AS4-003S**

• Stainless Steel

**AS4-004T**

• Titanium

#### Universal Chuck Handle



**AS4-005S**

• Stainless Steel

**AS4-005T**

• Titanium

## KNIVES and SPUDS

### KNIVES AND SPUDS

#### Tooke Corneal Knife



#### AS4-020

- Straight Blade With Rounded Cutting Edge

#### Gill Corneal Knife



#### AS4-021

- Slightly Curved 3.5mm Wide Blade

#### Hockey Epithelium Removal Knife



#### AS4-022

- Sharp 2.3mm Front Edge
- Semi-Sharp 4mm Bottom Edge

#### AS4-023

- Sharp 2.7 Front Edge
- Semi-Sharp 6mm Bottom Edge

#### Golf Club Foreign Body Spud



#### AS4-024

- Small Golf Club Shaped Tip With Sharp Edge
- For Foreign Body Removal

#### Francis Foreign Body Spud



#### AS4-025

- Flag Shaped Tip With Sharp Edge
- For Foreign Body Removal

## KNIVES and SPUDS

### Plange Foreign Body Spud



**AS4-026**

- Small Slightly Curved Spoon Shaped
- For Foreign Body Removal

### Double Ended Foreign Body Spud



**AS4-027**

- Double Ended Plange/Corneal Knife
- For Foreign Body Removal

## ALGER BRUSHES

### RUST RING REMOVAL INSTRUMENTS

### Alger Brush Rust Ring Remover



**AS4-100**

- Battery Powered Handpiece
- 0.5mm Burr

**AS4-101**

- Battery Powered Handpiece
- 1mm Burr

### Replacement Burr for Alger Brush



**AS4-100B**

- 0.5mm
- (Box of 5)

**AS4-101B**

- 1.0mm
- (Box of 5)

## PTERYGIUM INSTRUMENTS

### PTERYGIUM INSTRUMENTS

#### Round Diamond Burr with Power Handle



#### AS4-102

- Battery Powered Handpiece
- 2.5mm Diameter

#### AS4-103

- Battery Powered Handpiece
- 3.5mm Diameter

#### Round Diamond Burr (burr with chuck only)



#### AS4-102B

- 2.5mm Diameter



#### AS4-103B

- 3.5mm Diameter

#### Medium Grit Diamond with Power Handle



#### AS4-104

- 5.0mm Diameter

#### Medium Grit Diamond (burr with chuck only)



#### AS4-104B

- 5.0mm Diameter

## CORNEAL TREPHINES

### Corneal Trephine Blades



**AS4-300-4.0**  
• 4.0mm

**AS4-300-5.0**  
• 5.0mm

**AS4-300-6.0**  
• 6.0mm

**AS4-300-6.5**  
• 6.5mm

**AS4-300-7.0**  
• 7.0mm

**AS4-300-7.25**  
• 7.25mm

**AS4-300-7.5**  
• 7.5mm

**AS4-300-7.75**  
• 7.75mm

**AS4-300-8.0**  
• 8.0mm

**AS4-300-8.25**  
• 8.25mm

**AS4-300-8.5**  
• 8.5mm

**AS4-300-8.75**  
• 8.75mm

**AS4-300-9.0**  
• 9.0mm

**AS4-300-9.5**  
• 9.5mm

**AS4-300-13**  
• 13.0mm

**AS4-300-14**  
• 14.0mm