

ECO INSTRUMENTS

Speculums	50 – 1
Needle Holders.....	50 – 1-2
Scissors	50 – 2-4
Lens Loop	50 – 4
Capsulorrhexis Forceps.....	50 – 4-5
Utility Forceps	50 – 5
Corneal Forceps	50 – 6
Tying Forceps.....	50 – 7
Nucleus Rotator	50 – 7
Hooks.....	50 – 7-8
Lens Pusher.....	50 – 8
Choppers	50 – 9
Nucleus Manipulator	50 – 9

ECO INSTRUMENTS

SPECULUM

ECO Barraquer Wire Speculum

AS50-1001SR

- Adult Size
- 14mm Round Blades
- Standard Wire
- Stainless Steel

ECO Lieberman Temporal Speculum

AS50-1100KS

- Kratz Style Blades
- Adult 14mm
- Stainless Steel



AS50-1100VS

- V Style Blades
- Adult 14mm
- Stainless Steel

NEEDLE HOLDER

ECO Barraquer Needle Holder

AS50-11017S

- Standard Jaws 10mm
- Straight With Lock
- Stainless Steel



AS50-11018S

- Standard Jaws 10mm
- Straight Without Lock
- Stainless Steel

ECO INSTRUMENTS

ECO Barraquer Needle Holder



AS50-11019S

- Standard Jaws 10mm
- Curved With Lock
- Stainless Steel

AS50-11020S

- Standard Jaws 10mm
- Curved Without Lock
- Stainless Steel

SCISSORS

ECO Iris Scissors

AS50-12023S

- Pointed Tips
- Curved
- 105mm
- Stainless Steel

ECO Vannas Capsulotomy Scissors



AS50-12030S

- Straight
- Sharp Tips
- 6mm Blades
- Flat Handle
- Stainless Steel



AS50-12031S

- Curved
- Sharp Tips
- 6mm Blades
- Flat Handle
- Stainless Steel

ECO INSTRUMENTS

ECO Gills Scissors

AS50-12036

- Angled 8mm
- Sharp
- Stainless Steel

ECO Gills Capsulotomy Scissors

AS50-12039S

- Angled 7mm Blades
- Flat Handle
- Stainless Steel

ECO Westcott Tenotomy Scissors



AS50-12040S

- 13mm Blades
- Sharp
- Stainless Steel



AS50-12041S

- 16mm Blades
- Sharp
- Stainless Steel



AS50-12050S

- 13mm Blades
- Blunt
- Stainless Steel



AS50-12051S

- 16mm Blades
- Blunt
- Stainless Steel

ECO INSTRUMENTS

ECO Stevens Tenotomy Scissors

AS50-12055S

- Straight
- Blunt Tips
- Ring Handle
- Stainless Steel

AS50-12056S

- Curved
- Blunt Tips
- Ring Handle
- Stainless Steel

LENS LOOP

ECO Wilder Lens Loop with Serrations

AS50-14305S

- Stainless Steel

CAPSULORRHESIS FORCEPS

ECO Ultrata Capsulorrhesis Forceps

AS50-7100S

- Flat Handle
- 11.5mm Angled Shaft
- 85mm
- Stainless Steel

ECO INSTRUMENTS

ECO Ultrata Capsulorrhesis Forceps

AS50-7102S

- With Cystotome Tips
- Flat Handle
- 11.5mm Angled Shaft
- 85mm
- Stainless Steel



AS50-7103S

- With Cystotome Tips
- Long Flat Handle
- 11.5mm Angled Shaft
- Stainless Steel

AS50-7105S

- With Cystotome Tips
- Round Handle
- 11.5mm Angled Shaft
- 115mm
- Stainless Steel



AS50-7106S

- With Cystotome Tips
- Long Flat Handle
- 11.5mm Angled Shaft
- Stainless Steel

ECO Nugent Utility Forceps

AS50-7279S

- Serrated
- Stainless Steel

ECO INSTRUMENTS

CORNEAL FORCEPS

ECO Colibri Castroviejo Corneal Forceps

AS50-7301S

- .12mm
- 1 x 2 Teeth
- 85mm
- Stainless Steel

ECO Bonn Corneal Forceps

AS50-7320S

- Straight
- .12mm
- 1 x 2 Teeth
- 85mm
- Stainless Steel

ECO Castroviejo Corneal Forceps



AS50-7350S

- .12mm
- 1 x 2 Teeth
- 6mm Tying Platform
- Flat Handle
- Stainless Steel

AS50-7351S

- .3mm
- 1 x 2 Teeth
- 6mm Tying Platform
- Flat Handle
- Stainless Steel

AS50-7352S

- .5mm
- 1 x 2 Teeth
- 6mm Tying Platform
- Flat Handle
- Stainless Steel

ECO INSTRUMENTS

ECO McPherson Straight Tying Forceps



AS50-7430S

- 5mm Tying Platform
- Micro Handle 80mm
- Stainless Steel

ECO Kelman-McPherson Tying Forceps



AS50-7441S

- Angled 7mm
- Stainless Steel

AS50-7444S

- Angled 5mm
- Micro Handle 80mm
- Stainless Steel

ECO Bechert Nucleus Rotator

AS50-8019S

- Angled

ECO Kuglen Hooks



AS50-8022S

- Angled
- Stainless Steel

AS50-8023S

- Straight
- Stainless Steel

ECO INSTRUMENTS

ECO Lester Lens Pusher

AS50-8024S

- Angled
- Stainless Steel

AS50-8025S

- Straight
- Stainless Steel

ECO Sinskey Hook

**AS50-8030S**

- Angled
- .25mm Diameter
- Stainless Steel

ECO Jameson Muscle Hook

AS50-8050S

- .12mm
- Stainless Steel

AS50-8051S

- .10mm
- Stainless Steel

ECO INSTRUMENTS

CHOPPERS

ECO Nagahara I Phaco Chopper, RHD

**AS50-9020RS**

- To be Held in Left Hand
- 1mm Blunt Tip
- 0.5mm Cutting Edge
- Stainless Steel

ECOTriple Edged Chopper-Nichamin Style

**AS50-9031S**

- Multi Directional/Universal Cutting Edge
- Stainless Steel

ECO Quick Chopper with Rounded Edge-Seibel Type, RHD

**AS50-9036RS**

- To be Held in Left Hand
- Stainless Steel

ECO Drysdale Nucleus Manipulator

**AS50-9056S**

- Stainless Steel